

Product Datasheet

SSX2IP Antibody (orb1731323)

Description	Belongs to an adhesion system, which plays a role in the organization of homotypic, interneuronal an...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	A synthetic peptide of human SSX2IP
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 71 kDa; Observed MW: 71 kDa
Uniprot ID	Q9Y2D8
Entrez	117178
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.