

Product Datasheet

SGSM1 Antibody (orb1731365)

Description	Interacts with numerous Rab family members, functioning as Rab effector for some, and as GTPase acti...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	A synthetic peptide of human SGSM1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 130 kDa; Observed MW: 130 kDa
Uniprot ID	Q2NKQ1
Entrez	129049
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.