

Product Datasheet

Semaphorin 4D Antibody (orb1731375)

Description	Cell surface receptor for PLXN1B and PLXNB2 that plays an important role in cell-cell signaling. Pro...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	A synthetic peptide of human Semaphorin 4D
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 96 kDa; Observed MW: 120 kDa
Uniprot ID	Q92854
Entrez	10507
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.