

Product Datasheet

SCARB1 Antibody (orb1731385)

Description	Receptor for different ligands such as phospholipids, cholesterol ester, lipoproteins, phosphatidyls...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant protein of human Scavenging Receptor SR-BI
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 61 kDa; Observed MW: 80 kDa
Uniprot ID	Q8WTV0
Entrez	949
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.