

Product Datasheet

GTPase HRAS Antibody (orb1731432)

Description	Ras proteins bind GDP/GTP and possess intrinsic GTPase activity. Miscellaneous Mutations which change ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Recombinant protein of human RAS
Preservatives	50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Uniprot ID	P01111
Entrez	4893
Dilution Range	IP: 1/20; WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.