

Product Datasheet

RNA Polymerases I/II/III Subunit ABC1 Antibody (orb1731503)

Description	DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleosi...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human POLR2E
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Uniprot ID	P19388
Entrez	5434
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.