

Product Datasheet

PKC eta Antibody (orb1731513)

Description	PKC is activated by diacylglycerol which in turn phosphorylates a range of cellular proteins. PKC al...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human PKC eta
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 78 kDa; Observed MW: 78 kDa
Uniprot ID	P24723
Entrez	5583
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.