

Product Datasheet

SLC12A1 Antibody (orb1731613)

Description	Electrically silent transporter system. Mediates sodium and chloride reabsorption. Plays a vital rol...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human NKCC2
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 121 kDa; Observed MW: 150 kDa
Uniprot ID	Q13621
Entrez	6557
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.