

Product Datasheet

Neurofilament Heavy Polypeptide Antibody (orb1731624)

Description	Neurofilaments usually contain three intermediate filament proteins: L, M, and H which are involved ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Recombinant protein of human Neurofilament heavy polypeptide
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 112 kDa; Observed MW: 180-220 kDa
Uniprot ID	P12036
Entrez	4744
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.