

Product Datasheet

Myelin Protein Zero Antibody (orb1731661)

Description	Is an adhesion molecule necessary for normal myelination in the peripheral nervous system. It mediat...
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human Myelin Protein Zero
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Uniprot ID	P25189
Entrez	4359
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.