

Product Datasheet

LGALS3BP Antibody (orb1731745)

Description	Promotes integrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human LGALS3BP
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at - 20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 65 kDa; Observed MW: 100 kDa
Uniprot ID	Q08380
Entrez	3959
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.