

Product Datasheet

KDM1A Antibody (orb1731768)

Description	Histone demethylase that demethylates both 'Lys-4' (H3K4me) and 'Lys-9' (H3K9me) of histone H3, ther...
Species/Host	Rabbit
Reactivity	Hamster, Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	A synthetic peptide of human KDM1/LSD1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 93 kDa; Observed MW: 110 kDa
Uniprot ID	O60341
Entrez	23028
Dilution Range	IP: 1/20; WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.