

Product Datasheet

ICAM1 Antibody (orb1731825)

Description	ICAM1 a type I membrane protein of the immunoglobulin superfamily. Is a ligand for the leukocyte adh...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant protein of human ICAM1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 58 kDa; Observed MW: 89,92 kDa
Uniprot ID	P05362
Entrez	3383
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.