

Product Datasheet

Histone H2A Antibody (orb1731872)

Description	Core component of nucleosome which plays a central role in DNA double strand break (DSB) repair. Nuc...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of Histone H2A
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Uniprot ID	P04912
Entrez	852283
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.