

Product Datasheet

Calbindin Antibody (orb1732217)

Description	Buffers cytosolic calcium. May stimulate a membrane Ca ²⁺ -ATPase and a 3',5'-cyclic nucleotide phosph...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant protein of mouse Calbindin
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Uniprot ID	P12658
Entrez	12307
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.