

Product Datasheet

AP3M1 Antibody (orb1732307)

Description	Part of the AP-3 complex, an adaptor-related complex which is not clathrin-associated. The complex i...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human AP3M1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Uniprot ID	Q9Y2T2
Entrez	26985
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.