

Product Datasheet

Glutamine Synthetase Antibody (orb1732447)

Description	The protein encoded by this gene belongs to the glutamine synthetase family. It catalyzes the synthe...
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human Glutamine Synthetase
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Uniprot ID	P15104
Entrez	2752
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.