

Product Datasheet

CKMT1 Antibody (orb1732562)

Description	Creatine kinase MT is responsible for the transfer of high energy phosphate from mitochondria to the...
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human CKMT1 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Uniprot ID	P12532
Entrez	1159
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.