

Product Datasheet

Galectin 3 Antibody (orb1732582)

Description	Galactose-specific lectin which binds IgE. May mediate with the alpha-3, beta-1 integrin the stimula...
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic Peptide of Galectin-3
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Uniprot ID	P17931
Entrez	3958
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.