

Product Datasheet

PGP9.5 Antibody (orb1732630)

Description	PGP9.5 is a single chain transmembraneous glycoprotein which represents the low molecular weight iso...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Purified recombinant human UCHL1 / PGP9.5 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Uniprot ID	P09936
Entrez	7345
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.