

Product Datasheet

ICAM1 Antibody (orb1732647)

Description	ICAM1 a type I membrane protein of the immunoglobulin superfamily. Is a ligand for the leukocyte adh...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Purified recombinant human CD54(ICAM-1) protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	Calculated MW: 58 kDa; Observed MW: 96 kDa
Uniprot ID	P05362
Entrez	3383
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.