

Product Datasheet

Transcription Initiation Factor IIB Antibody (orb1732654)

Description	General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA poly...
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human TFIIIB protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Uniprot ID	Q00403
Entrez	2959
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.