

Product Datasheet

LC3A/B Antibody (orb1732704)

Description	Macroautophagy is the major inducible pathway for the general turnover of cytoplasmic constituents i...
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to human LC3B protein
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	Calculated MW: 14 kDa; Observed MW: 14,16 kDa
Uniprot ID	Q9H492, Q9GZQ8
Entrez	84557, 81631
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.