

Product Datasheet

Toll-Like Receptor 3 Antibody (orb1732714)

Description	The protein encoded by this gene is a member of the Toll-like receptor (TLR) family which plays a fu...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human Toll-Like Receptor 3 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	Calculated MW: 104 kDa; Observed MW: Refer to figures
Uniprot ID	O15455
Entrez	7098
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.