

Product Datasheet

Wnt3 Antibody (orb1732718)

Description	Ligand for members of the frizzled family of seven transmembrane receptors. Wnt-3 and Wnt-3a play di...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human Wnt3 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Uniprot ID	P56703
Entrez	7473
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.