

Product Datasheet

NUP98 Antibody (orb1732795)

Description	Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. Nup98 and Nup96 are invo...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Purified recombinant human NUP98 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 198 kDa; Observed MW: 98 kDa
Uniprot ID	P52948
Entrez	4928
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.