

Product Datasheet

CDK5 Antibody (orb1732810)

Description	Activated by cyclins but by p35 (CDK5R1) and p39. An important regulator of neuronal positioning dur...
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Purified recombinant human CDK5(N-terminus) protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 33 kDa; Observed MW: 36 kDa
Uniprot ID	Q00535
Entrez	1020
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.