

Product Datasheet

USP7 Antibody (orb1732827)

Description	Hydrolase that deubiquitinates target proteins such as FOXO4, p53/TP53, MDM2, ERCC6, DNMT1, UHRF1, P...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human HAUSP / USP7 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 128 kDa; Observed MW: 135 kDa
Uniprot ID	Q93009
Entrez	7874
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.