

Product Datasheet

CD68 Antibody (orb1732859)

Description	CD68 belongs to a family of acidic, highly glycosylated lysosomal glycoproteins (LGPs) that includes...
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC
Immunogen	Synthetic Peptide of CD68
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P34810
Entrez	968
Dilution Range	IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.