

Product Datasheet

COL1A1 Antibody (orb1732884)

Description	<span 46="" 479="" 907="" 959"="" data-label="Page-Footer" style="font-family:arial, helvetica, clean, sans-serif;font-size:13px;background-color:#FFFFFFF...</td></tr><tr><td>Species/Host</td><td>Rabbit</td></tr><tr><td>Reactivity</td><td>Human</td></tr><tr><td>Conjugation</td><td>Unconjugated</td></tr><tr><td>Tested Applications</td><td>ICC, WB</td></tr><tr><td>Immunogen</td><td>Recombinant protein of human Collagen I</td></tr><tr><td>Preservatives</td><td>50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin</td></tr><tr><td>Concentration</td><td>1mg/ml</td></tr><tr><td>Storage</td><td>Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.</td></tr><tr><td>Note</td><td>For research use only</td></tr><tr><td>Isotype</td><td>Rabbit IgG</td></tr><tr><td>Clonality</td><td>Monoclonal</td></tr><tr><td>MW</td><td>139 kDa</td></tr><tr><td>Entrez</td><td>1277</td></tr><tr><td>Dilution Range</td><td>WB: 1/1000; ICC: 1/50</td></tr><tr><td>Expiration Date</td><td>12 months from date of receipt.</td></tr></table></div><div data-bbox="> Biorbyt Ltd. 7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom Email: info@biorbyt.com Phone: +44 (0) 1223 859-353 Fax: +44 (0)1223 280 240
--------------------	---