

Product Datasheet

LC3A Antibody (orb1732907)

Description	Autophagy marker Light Chain 3 (LC3) was originally identified as a subunit of microtubule-associate...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Purified recombinant protein expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 14 kDa; Observed MW: 14,16 kDa
Uniprot ID	Q9H492
Entrez	84557
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.