

Product Datasheet

Cystatin C Antibody (orb1732910)

Description	Cystatin C is a 14 kDa member of the Cystatin superfamily of cysteine protease inhibitors. Most cell...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant protein expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Uniprot ID	P01034
Entrez	1471
Dilution Range	WB: 1/500-1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.