

Product Datasheet

HspB8 Antibody (orb1732911)

Description	HSPB8 (HSP22) is a member of the small heat shock protein superfamily and the human protein is most ...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant protein expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Uniprot ID	Q9UJY1
Entrez	26353
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.