

Product Datasheet

beta I Tubulin Antibody (orb1732914)

Description	Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable ...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic Peptide of β I tubulin
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Uniprot ID	Q9H4B7
Entrez	81027
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.