

Product Datasheet

gamma Tubulin Antibody (orb1732956)

Description	TUBG1 Tubulin is the major constituent of microtubules. Gamma tubulin is found at microtubule organi...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide conjugated to KLH.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P23258
Entrez	7283
Dilution Range	IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.