

Product Datasheet

IL-8 Antibody (orb1732964)

Description	IL-8 is a chemotactic factor that attracts neutrophils, basophils, and T-cells, but not monocytes. I...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide conjugated to KLH.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P10145
Entrez	3576
Dilution Range	IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.