

Product Datasheet

Phospho-TAK1 (Ser439) Antibody (orb1732988)

Description	Component of a protein kinase signal transduction cascade. Mediator of TRAF6 and TGF-beta signal tra...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser439 of human TAK1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 67 kDa; Observed MW: 78 kDa
Uniprot ID	043318
Entrez	6885
Dilution Range	IP: 1/20; WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.