

Product Datasheet

Phospho-Histone H2A.X (Ser139) Antibody (orb1732991)

Description

Variant histone H2A which replaces

nucleosomes. Nucleosomes wrap and...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IP, WB

Applications
Immunogen

A synthetic phosphopeptide corresponding to residues surrounding Ser139 of human gamma H2A.X

Preservatives

50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin

Concentration

1mg/ml

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Monoclonal

MW

Calculated MW: 15 kDa; Observed MW: 15 kDa

Uniprot ID

P16104

Entrez

3014

Dilution Range

IP: 1/20; WB: 1/500-1/1000

Expiration Date

12 months from date of receipt.