

Product Datasheet

Phospho-Smad3 (Ser423/425) Antibody (orb1732995)

Description	Smad3 transcription factor phosphorylated and activated by TGF-beta-type receptors. A receptor-regul...
Species/Host	Rabbit
Reactivity	Hamster, Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser423/Ser425 of human Smad3
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 48 kDa; Observed MW: 52 kDa
Uniprot ID	P84022
Entrez	4088
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.