

Product Datasheet

Phospho-Tyrosine Antibody (orb1733001)

Description	Tyrosine phosphorylation plays a key role in cellular signaling. Research studies have shown that in...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Purified Protein
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Refer to figures
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.