

Product Datasheet

Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) Antibody (orb1733004)

Description	Serine/threonine protein kinase that plays a role in a variety of different signaling pathways inclu...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602)
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 64,80,75 kDa; Observed MW: 64,80,75 kDa
Uniprot ID	Q9P286, O96013, Q9NQU5
Entrez	57144, 56924, 10298
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.