

Product Datasheet

Phospho-TrkB (Tyr705) Antibody (orb1733005)

Description	The family of Trk receptor tyrosine kinases consists of TrkA, TrkB and TrkC. While the sequence of t...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Tyr705 of human TrkB
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 92 kDa; Observed MW: 140,90 kDa
Uniprot ID	Q16620
Entrez	4915
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.