

Product Datasheet

Phospho-Src (Tyr530) Antibody (orb1733007)

Description	Src activity is regulated by tyrosine phosphorylation at two sites, but with opposing effects. While...
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Tyr529 of human Src
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Uniprot ID	P12931
Entrez	6714
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.