

Product Datasheet

Phospho-RSK1 p90 (Ser380) Antibody (orb1733010)

Description	Rsk1 is a member of a family of 90kDa ribosomal protein S6 kinases, which includes Rsk1, Rsk2 and Rs...
Species/Host	Rabbit
Reactivity	Hamster, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser380 of human RSK1 p90
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 83 kDa; Observed MW: 90 kDa
Uniprot ID	Q15418
Entrez	6195
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.