

Product Datasheet

Phospho-PKC delta (Ser299) Antibody (orb1733019)

Description

Protein kinase C (PKC) is a family of

kinases that can be act...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, IP, WB

Applications

Immunogen

A synthetic phosphopeptide corresponding to residues surrounding Ser299 of human PKC delta

Preservatives

50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin

Concentration

1mg/ml

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Monoclonal

MW

Calculated MW: 78 kDa; Observed MW: 78 kDa

Uniprot ID

[Q05655](#)

Entrez

[5580](#)

Dilution Range

IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100

Expiration Date

12 months from date of receipt.