

Product Datasheet

Phospho-FAK (Tyr576) Antibody (orb1733036)

Description	This gene encodes a cytoplasmic protein tyrosine kinase which is found concentrated in the focal adh...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Tyr576 of human FAK
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 119 kDa; Observed MW: 125 kDa
Uniprot ID	Q05397
Entrez	5747
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.