

Product Datasheet

Phospho-delta 1 Catenin (Ser268) Antibody (orb1733044)

Description

Catenin δ -1 (p120 catenin) has an

followed by a regulatory domain...

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic phosphopeptide corresponding to residues surrounding Ser268 of human δ 1 Catenin/CAS

Preservatives

50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin

Concentration

1mg/ml

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Monoclonal

MW

Calculated MW: 108 kDa; Observed MW: 100 kDa

Uniprot ID

O60716

Entrez

1500

Dilution Range

WB: 1/500-1/1000

Expiration Date

12 months from date of receipt.