

Product Datasheet

Phospho-PRKD2 (Ser876) Antibody (orb1733071)

Description	Converts transient diacylglycerol (DAG) signals into prolonged physiological effects, downstream of ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser876 of human PKD2
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 97 kDa; Observed MW: 105 kDa
Uniprot ID	Q9BZL6
Entrez	25865
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.