

Product Datasheet

Phospho-MEK1 (Thr292) Antibody (orb1733083)

Description	The protein encoded by this gene is a member of the dual specificity protein kinase family, which ac...
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Thr292 of human MEK1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Uniprot ID	Q02750
Entrez	5604
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.