

Product Datasheet

Phospho-EGFR (Tyr1068) Antibody (orb1733097)

Description	EGFR is a receptor tyrosine kinase. Receptor for epidermal growth factor (EGF) and related growth fa...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Tyr1068 of human EGFR
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 134 kDa; Observed MW: 175 kDa
Uniprot ID	P00533
Entrez	1956
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.